

DOC-001

DRUG MASTER FILE (DMF)

MODULE 3.2.S: DRUG SUBSTANCE (API)

Version: 1.0

Date: [To be determined]

Author: Stefan Schönwälder

This document is for demonstration and training purposes only.

MODULE 3.2.S: DRUG SUBSTANCE (API)

Fampridin (4-Aminopyridine)

Submission to: Pharmaceuticals and Medical Devices Agency (PMDA), Japan

DMF Type: Type II Drug Master File

Format: Common Technical Document (CTD)

DMF Number: DMF-FAM-001

Submission Date: [To be determined]

Applicant: [Manufacturing Site Name]

Address: [Site Address]

Document Control:

Version	Date	Author	Description
1.0	04.12.2025	Regulatory Affairs	Initial DMF for PMDA submission

TABLE OF CONTENTS

MODULE 3.2.S: DRUG SUBSTANCE

- 3.2.S.1 General Information
 - 3.2.S.1.1 Nomenclature
 - 3.2.S.1.2 Structure

- 3.2.S.1.3 General Properties

- 3.2.S.2 Manufacture
 - 3.2.S.2.1 Manufacturer(s)
 - 3.2.S.2.2 Description of Manufacturing Process and Process Controls
 - 3.2.S.2.3 Control of Materials
 - 3.2.S.2.4 Controls of Critical Steps and Intermediates
 - 3.2.S.2.5 Process Validation and/or Evaluation
 - 3.2.S.2.6 Manufacturing Process Development

- 3.2.S.3 Characterization
 - 3.2.S.3.1 Elucidation of Structure and Other Characteristics
 - 3.2.S.3.2 Impurities

- 3.2.S.4 Control of Drug Substance
 - 3.2.S.4.1 Specification
 - 3.2.S.4.2 Analytical Procedures
 - 3.2.S.4.3 Validation of Analytical Procedures
 - 3.2.S.4.4 Batch Analyses
 - 3.2.S.4.5 Justification of Specification

- 3.2.S.5 Reference Standards or Materials

- 3.2.S.6 Container Closure System

- 3.2.S.7 Stability
 - 3.2.S.7.1 Stability Summary and Conclusions

- 3.2.S.7.2 Post-approval Stability Protocol and Stability Commitment
- 3.2.S.7.3 Stability Data

3.2.S.1 GENERAL INFORMATION

3.2.S.1.1 NOMENCLATURE

International Nonproprietary Name (INN): Fampridine (WHO)

Alternative INN: Fampridin (German, European)

Chemical Name (IUPAC): 4-Aminopyridine

CAS Registry Number: 504-24-5

Molecular Formula: C₅H₆N₂

Molecular Weight: 94.11 g/mol

Compendial Status:


- **DAC (Deutscher Arzneimittel-Codex):** Monograph available
- **Japanese Pharmacopoeia (JP):** Not currently listed (to be established)

Pharmacological Class: Potassium Channel Blocker

Therapeutic Use: Treatment of walking impairment in patients with Multiple Sclerosis

3.2.S.1.2 STRUCTURE

Structural Formula:



Chemical Structure (2D):

- Pyridine ring with amino group at position 4

Stereochemistry: Not applicable (no chiral centers)

Molecular Formula: C₅H₆N₂

Structural Confirmation:

The structure of Fampridine has been confirmed by the following analytical techniques:

- **Infrared Spectroscopy (IR):** Characteristic absorption bands consistent with amino group and aromatic ring
- **Nuclear Magnetic Resonance (¹H-NMR and ¹³C-NMR):** Confirmed substitution pattern on pyridine ring
- **Mass Spectrometry (MS):** Molecular ion peak at m/z = 94.11 (consistent with molecular formula)
- **Elemental Analysis:** C, H, N content consistent with C₅H₆N₂

3.2.S.1.3 GENERAL PROPERTIES

Physical Form: White to off-white crystalline powder

Solubility:

- Water: Freely soluble (>100 mg/mL at 25°C)
- Ethanol: Soluble (~50 mg/mL at 25°C)
- Methanol: Soluble
- Acetone: Sparingly soluble
- Dichloromethane: Practically insoluble

pKa: 9.2 (basic amine group)

Melting Point: 158-160°C (literature value)

Partition Coefficient (log P): -0.20 (octanol/water at pH 7.4) – Hydrophilic compound

Hygroscopicity: Non-hygroscopic (water uptake <0.2% at 25°C/80% RH)

Polymorphism: No polymorphs identified. Single crystalline form (Form I) consistently obtained from manufacturing process.

3.2.S.2 MANUFACTURE

3.2.S.2.1 MANUFACTURER(S)

API Manufacturer:

- **Name:** [Manufacturing Site Name]
- **Address:** [Site Address, City, Country]

- **GMP Certification:** EU GMP Certificate issued by [Authority], Certificate No. [XXX], valid until [Date]
- **Manufacturing Authorization:** [Authorization Number]

Regulatory Status:

- The manufacturing site is inspected and approved for API production under EU GMP guidelines.
- PMDA pre-approval inspection requested (pending).

3.2.S.2.2 DESCRIPTION OF MANUFACTURING PROCESS AND PROCESS CONTROLS

3.2.S.2.2.1 MANUFACTURING PROCESS OVERVIEW

Process Type: Chemical Synthesis

Number of Steps: 6 main unit operations

Batch Size: 100 kg (commercial scale)

Process Flow Diagram:

```
Raw Materials
  ↓
[Step 1: REACTION] → Synthesis Reactor R-101
  ↓
[Step 2: WORKUP] → Neutralization, Extraction
  ↓
[Step 3: CRYSTALLIZATION] → Purification
  ↓
[Step 4: FILTRATION] → Solid/Liquid Separation
  ↓
[Step 5: DRYING] → Vacuum Dryer VD-101
  ↓
```

[Step 6: MILLING] → Particle Size Control

↓

[PACKAGING] → Polyethylene-lined Drums

↓

API (Fampridin)

3.2.S.2.2.2 DETAILED MANUFACTURING PROCESS

STEP 1: REACTION (SYNTHESIS)

Equipment: Glass-lined reactor R-101 (500 L capacity)

Starting Materials:

- Starting Material A: [Chemical Name], CAS [XXX-XX-X]
- Starting Material B: [Chemical Name], CAS [XXX-XX-X]
- Solvent: Ethanol (Ph. Eur. grade)

Procedure:

1. Charge Starting Material A ([XX] kg) and Starting Material B ([XX] kg) to reactor
2. Add ethanol ([XX] L)
3. Heat to **80 ± 2°C** under nitrogen atmosphere
4. Maintain temperature for **4.0 ± 0.5 hours** with agitation
5. Monitor pH: Target **7.0 ± 0.5**

Critical Process Parameters (CPPs):

- **Reaction Temperature:** 80 ± 2°C (monitored continuously via calibrated temperature probe)
- **Reaction Time:** 4.0 ± 0.5 hours

- **pH:** 7.0 ± 0.5

In-Process Control (IPC):

- **IPC-1:** HPLC purity at end of reaction: $\geq 95\%$
- If purity $< 95\%$: Extend reaction time (maximum 5 hours total)

Theoretical Yield (Step 1): $\sim 95\%$

STEP 2: WORKUP (NEUTRALIZATION AND EXTRACTION)

Equipment: Workup reactor R-102 (300 L), Separator vessel

Procedure:

1. Transfer reaction mixture to workup reactor R-102
2. Cool to $25 \pm 5^\circ\text{C}$
3. Adjust pH to **6.5-7.5** using [Acid/Base]
4. Add extraction solvent ([XX] L)
5. Mix for 30 minutes, allow phase separation
6. Separate aqueous and organic layers
7. Collect organic layer containing Fampridin

Critical Process Parameters:

- **pH after neutralization:** 6.5-7.5 (monitored via calibrated pH meter)

In-Process Control:

- **IPC-2:** pH verification: 6.5-7.5
- Action if OOS: Adjust with acid or base

STEP 3: CRYSTALLIZATION (PURIFICATION)

Equipment: Crystallizer C-101 (jacketed vessel, 400 L)

Procedure:

1. Transfer organic phase to crystallizer
2. Concentrate by distillation under vacuum (remove ~50% solvent)
3. Cool to **$10 \pm 5^{\circ}\text{C}$** to induce crystallization
4. Hold at crystallization temperature for **2 ± 0.5 hours**
5. Age crystals with agitation

Critical Process Parameters:

- **Crystallization Temperature:** $10 \pm 5^{\circ}\text{C}$
- **Crystallization Hold Time:** 2 ± 0.5 hours

In-Process Control:

- **IPC-3:** HPLC purity after crystallization: $\geq 98\%$
- If purity $< 98\%$: Recrystallize

Expected Yield (Step 3): ~90% (from Step 1 input)

STEP 4: FILTRATION

Equipment: Nutsche filter F-101 (filter area: 1 m²)

Procedure:

1. Filter crystals under vacuum
2. Wash with cold ethanol ([XX] L)
3. Dry filter cake under vacuum for 30 minutes

In-Process Control: Visual inspection (crystals white to off-white)

STEP 5: DRYING

Equipment: Vacuum dryer VD-101 (tray dryer, 200 L capacity)

Procedure:

1. Transfer wet filter cake to vacuum dryer
2. Dry at **60 ± 5°C** under vacuum **<50 mbar**
3. Dry until Loss on Drying (LOD) **≤0.50%**
4. Typical drying time: 8-12 hours
5. Cool to room temperature under nitrogen

Critical Process Parameters:

- **Drying Temperature:** 60 ± 5°C

- **Vacuum Level:** <50 mbar

In-Process Control:

- **IPC-4:** Loss on Drying (LOD) $\leq 0.50\%$

- Action if >0.50%: Continue drying

Expected Yield (Step 5): ~88% (from Step 1 input)

STEP 6: MILLING

Equipment: Impact mill M-101

Procedure:

1. Mill dried API to achieve target particle size
2. Pass through sieve (if applicable)

Critical Process Parameter:

- **Target Particle Size Distribution:** D50: 20-80 μm

In-Process Control:

- **IPC-5:** Particle size distribution by laser diffraction: D50: 20-80 μm

- If out of range: Re-mill or blend with appropriate batch

Final Yield (Overall): 85-95% (from Step 1 starting materials)

STEP 7: PACKAGING

Equipment: Packaging area (ISO 8 cleanroom)

Procedure:

1. Transfer milled API to polyethylene-lined drums (25 kg per drum)
2. Close drums with secure lids and tamper-evident seals
3. Label drums with:
 - Product name: Fampridin
 - Batch number
 - Manufacturing date
 - Quantity
 - Storage conditions: "Store at 15-25°C, protect from light and moisture"

Packaging Material: Food-grade polyethylene liner, steel drum (exterior)

3.2.S.2.3 CONTROL OF MATERIALS

3.2.S.2.3.1 RAW MATERIALS

All raw materials are purchased from qualified suppliers and tested per specifications before use.

Starting Material A:

- **Specification:** [Details of identity, purity, impurity limits]
- **Supplier:** [Supplier Name, Country]
- **Testing:** Identity (IR), Assay (HPLC), Impurities (HPLC), Water Content (KF)

Starting Material B:

- **Specification:** [Details]
- **Supplier:** [Supplier Name, Country]
- **Testing:** Identity, Assay, Impurities

Solvents (Ethanol):

- **Specification:** Ph. Eur. grade, $\geq 99.5\%$
- **Supplier:** [Supplier Name]
- **Testing:** Identity (GC), Assay (GC), Water Content (KF)

Supplier Qualification:

- All suppliers are qualified per SOP-SUPPLIER-001
- Supplier audits conducted (on-site or questionnaire)
- Quality agreements in place

3.2.S.2.3.2 EXCIPIENTS / PROCESSING AIDS

Not Applicable: Fampridin is an API. No excipients are added during synthesis.

Processing Aids (Solvents) Removed:

- Ethanol is removed during crystallization and drying steps

- Residual solvent content controlled per specification (≤ 5000 ppm, per ICH Q3C Class 3 solvent)

3.2.S.2.4 CONTROLS OF CRITICAL STEPS AND INTERMEDIATES

Critical Steps are process steps where parameters directly impact Critical Quality Attributes (CQAs).

Summary of Critical Steps:

Step	Critical Parameter	Target/Range	Linked CQA	Control Method
Reaction	Temperature	$80 \pm 2^{\circ}\text{C}$	Assay, Impurities	Continuous monitoring (calibrated probe)
Reaction	Time	4.0 ± 0.5 h	Assay, Impurities	Timer
Reaction	pH	7.0 ± 0.5	Assay, Impurities	Continuous pH probe
Workup	pH	6.5-7.5	Impurities	pH meter
Crystallization	Temperature	$10 \pm 5^{\circ}\text{C}$	Particle Size, Purity	Temperature probe
Crystallization	Time	2 ± 0.5 h	Particle Size, Purity	Timer
Drying	Temperature	$60 \pm 5^{\circ}\text{C}$	Water Content, Residual Solvents	Temperature probe
Drying	Vacuum	<50 mbar	Water Content, Residual Solvents	Vacuum gauge
Milling	Particle Size	D50: 20-80 μm	Particle Size, Dissolution	Laser Diffraction

All critical parameters are monitored and recorded in the Batch Manufacturing Record.

Intermediates:

- No isolated intermediates are stored or transferred between sites
- All synthesis steps occur at the same manufacturing site in a continuous process

3.2.S.2.5 PROCESS VALIDATION AND/OR EVALUATION

Validation Status: Process is validated per ICH Q7 and EU GMP Annex 15

Validation Approach: Prospective validation using **3 consecutive commercial-scale batches**

Validation Protocol: PVP-FAM-001 (Process Validation Protocol – Fampridin)

Validation Batches:

- Batch 1: [Batch Number]
- Batch 2: [Batch Number]
- Batch 3: [Batch Number]

Validation Results Summary:

Parameter	Batch 1	Batch 2	Batch 3	Mean	RSD (%)	Acceptance Criteria	Pass/Fail
Assay (%)	99.8	99.9	99.7	99.8	0.10	99.0-101.0	Pass
Total Impurities (%)	0.32	0.29	0.35	0.32	9.4	≤0.50	Pass
Water Content (%)	0.25	0.30	0.28	0.28	9.0	≤0.50	Pass
Particle Size D50 (µm)	45	50	48	48	5.2	20-80	Pass
Yield (%)	90	88	91	90	1.7	85-95	Pass

Conclusion: The manufacturing process consistently produces Fampridin API meeting all specifications. Process is validated.

Validation Report: PVR-FAM-001 (on file)

3.2.S.2.6 MANUFACTURING PROCESS DEVELOPMENT

Development History:

Phase 1: Laboratory Scale (2022-2023)

- Synthesis route established
- Critical Quality Attributes (CQAs) identified
- Critical Process Parameters (CPPs) determined through Design of Experiments (DoE)

Phase 2: Pilot Scale (2023-2024)

- Scale-up from 1 kg to 10 kg batches
- Process optimization (crystallization conditions, drying parameters)
- Impurity profile established

Phase 3: Commercial Scale (2024-2025)

- Scale-up to 100 kg batches
- Process validation (3 batches)
- Qualification of manufacturing equipment (IQ/OQ/PQ)

Process Changes During Development:

- **Change 1 (2023):** Crystallization solvent modified (from methanol to ethanol) to reduce residual solvent limits

- **Change 2 (2024):** Drying temperature reduced from 70°C to 60°C to minimize thermal degradation

Current Process: Locked and validated (no further changes planned)

3.2.S.3 CHARACTERIZATION

3.2.S.3.1 ELUCIDATION OF STRUCTURE AND OTHER CHARACTERISTICS

Structural Elucidation:

The structure of Fampridin has been confirmed by comprehensive analytical characterization:

1. Infrared Spectroscopy (IR):

- Characteristic absorption bands:
 - N-H stretch: 3300-3400 cm^{-1} (primary amine)
 - C=C aromatic stretch: 1580-1600 cm^{-1}
 - C-N stretch: 1250-1350 cm^{-1}
- IR spectrum matches Fampridin reference standard

2. Nuclear Magnetic Resonance (NMR):

- **^1H -NMR (DMSO- d_6):** Consistent with 4-aminopyridine structure
 - Aromatic protons (3H): δ 7.5-8.0 ppm
 - Amino protons (2H): δ 5.8 ppm (broad)
- **^{13}C -NMR:** Five carbon signals consistent with substituted pyridine

3. Mass Spectrometry (MS):

- **ESI-MS:** $[M+H]^+$ = 95.1 (molecular ion peak)
- Fragmentation pattern consistent with Fampridin

4. Elemental Analysis:

- Theoretical ($C_5H_6N_2$): C = 63.81%, H = 6.43%, N = 29.76%
- Found: C = 63.75%, H = 6.50%, N = 29.70%
- Within $\pm 0.4\%$ of theoretical values

Polymorphism:

- Powder X-Ray Diffraction (PXRD) analysis performed on 10 batches
- Single crystalline form (Form I) consistently obtained
- No polymorphs identified

Particle Size and Morphology:

- Scanning Electron Microscopy (SEM): Crystalline particles, irregular shape
- Typical particle size (D50): 40-50 μm

3.2.S.3.2 IMPURITIES

3.2.S.3.2.1 POTENTIAL IMPURITIES

Impurities have been identified through:

- Process understanding (reaction mechanism, degradation pathways)
- Forced degradation studies (acid, base, oxidative, thermal, photolytic stress)

- Analysis of manufacturing batches

Types of Impurities:

1. Organic Impurities (Process-Related):

- **Impurity A:** [Chemical Name/Structure] – Starting material residue
- **Impurity B:** [Chemical Name/Structure] – Reaction by-product
- **Impurity C:** [Chemical Name/Structure] – Degradation product

2. Inorganic Impurities:

- Residual metals (from reactors): Controlled by equipment materials (stainless steel, glass-lined)
- Specification: Heavy metals ≤ 10 ppm (if tested)

3. Residual Solvents:

- **Ethanol:** Class 3 solvent per ICH Q3C
- Specification: ≤ 5000 ppm
- Typical levels: < 500 ppm (well below limit)

3.2.S.3.2.2 IMPURITY QUALIFICATION

Impurity Qualification Thresholds (per ICH Q3A(R2)):

Daily Dose	Reporting Threshold	Identification Threshold	Qualification Threshold
≤ 2 g/day	0.05%	0.10% or 1.0 mg/day	0.15% or 1.0 mg/day

Fampridin Typical Dose: 20 mg/day (for oral formulation)

Since dose <2 g/day:

- Reporting: $\geq 0.05\%$
- Identification: $\geq 0.10\%$ or ≥ 1.0 mg intake/day
- Qualification: $\geq 0.15\%$ or ≥ 1.0 mg intake/day

Impurity Profile in Commercial Batches:

Impurity	Typical Level (%)	Maximum Observed (%)	Identification Required?	Qualification Status
Impurity A	0.05	0.08	No (<0.10%)	Not Required
Impurity B	0.10	0.12	Yes ($\geq 0.10\%$)	Qualified (toxicological assessment on file)
Impurity C	0.03	0.05	No (<0.10%)	Not Required
Total	0.18	0.25	N/A	N/A

All impurities are controlled within specification limits (Any single $\leq 0.10\%$, Total $\leq 0.50\%$).

3.2.S.4 CONTROL OF DRUG SUBSTANCE**3.2.S.4.1 SPECIFICATION****Fampridin API Specification (SPEC-FAM-001):**

Test	Method	Acceptance Criteria
Appearance	Visual	White to off-white crystalline powder
Identity (IR)	Ph. Eur. 2.2.24	Conforms to reference standard
Identity (HPLC)	In-house HPLC-001	Retention time matches

		reference $\pm 2\%$
Assay (Content)	HPLC-001	99.0 - 101.0% (on anhydrous basis)
Impurities	HPLC-001	Any single impurity: $\leq 0.10\%$ Total impurities: $\leq 0.50\%$
Water Content	Ph. Eur. 2.5.12 (Karl Fischer)	$\leq 0.50\%$
Residual Solvents (Ethanol)	GC-001	≤ 5000 ppm
Particle Size Distribution	Laser Diffraction	D50: 20-80 μm
Microbial Limits	Ph. Eur. 2.6.12	Total aerobic count: ≤ 1000 CFU/g Yeast/Mold: ≤ 100 CFU/g E. coli: Absent Salmonella: Absent

All tests are performed on each batch before release.

3.2.S.4.2 ANALYTICAL PROCEDURES

Summary of Analytical Methods:

Test	Method Reference	Method Type
Appearance	Visual inspection	Pharmacopoeial (general)
Identity (IR)	IR Spectroscopy	Pharmacopoeial (Ph. Eur. 2.2.24)
Identity (HPLC)	HPLC-001	In-house (non-compendial)
Assay	HPLC-001	In-house (non-compendial)
Impurities	HPLC-001	In-house (non-compendial)
Water Content	Karl Fischer Titration	Pharmacopoeial (Ph. Eur. 2.5.12)
Residual Solvents	GC-001	In-house (based on USP <467>)
Particle Size	Laser Diffraction	In-house (ISO 13320)
Microbial Limits	Ph. Eur. 2.6.12	Pharmacopoeial

Detailed analytical procedures are provided in Appendix 1 (attached separately).

3.2.S.4.3 VALIDATION OF ANALYTICAL PROCEDURES

All non-compendial analytical methods are validated per ICH Q2(R1).

Validation Summary:

HPLC Method (Assay and Impurities):

- **Validation Protocol:** AMV-HPLC-001
- **Validation Report:** AMVR-HPLC-001 (on file)
- **Parameters Validated:**
 - Specificity: ✓ (forced degradation studies, resolution ≥ 2.0)
 - Linearity: ✓ ($R^2 = 0.9998$, range 50-150%)
 - Accuracy: ✓ (Recovery 99.5-100.8%, RSD <2%)
 - Precision: ✓ (Repeatability RSD 0.8%, Intermediate Precision RSD 1.2%)
 - LOD: 0.03%, LOQ: 0.05%
 - Robustness: ✓ (method stable to minor variations)

GC Method (Residual Solvents):

- **Validation Protocol:** AMV-GC-001
- **Validation Report:** AMVR-GC-001 (on file)
- **Parameters Validated:** Specificity, Linearity ($R^2 > 0.999$), Accuracy (Recovery 95-105%), Precision (RSD <5%), LOD/LOQ

All analytical methods are validated and suitable for their intended purpose.

3.2.S.4.4 BATCH ANALYSES

Batch Analysis Data (Representative Commercial Batches):

Batch Number	Manufacture Date	Assay (%)	Total Impurities (%)	Water (%)	Ethanol (ppm)	D50 (µm)	Microbial (CFU/g)
2025-01-001	Jan 2025	99.8	0.32	0.25	450	45	<10
2025-01-002	Jan 2025	99.9	0.29	0.30	380	50	<10
2025-02-001	Feb 2025	99.7	0.35	0.28	520	48	<10
2025-03-001	Mar 2025	99.9	0.28	0.32	410	52	<10
2025-04-001	Apr 2025	99.8	0.30	0.26	460	46	<10
2025-05-001	May 2025	99.8	0.33	0.29	490	49	<10

Conclusion: All batches meet specification. Process is consistent and under control.

Certificates of Analysis (CoA) for all batches are available upon request.

3.2.S.4.5 JUSTIFICATION OF SPECIFICATION

Specification limits are justified based on:

1. Regulatory Guidelines:

- ICH Q6A: Specifications – Test Procedures and Acceptance Criteria for New Drug Substances
- ICH Q3A(R2): Impurities in New Drug Substances

- ICH Q3C(R8): Impurities – Guideline for Residual Solvents

2. Manufacturing Data:

- Batch analysis data from development and commercial batches (n=30 batches)
- All batches consistently meet proposed specification limits

3. Stability Data:

- Long-term stability studies (25°C/60% RH, up to 24 months)
- Accelerated stability studies (40°C/75% RH, 6 months)
- No significant degradation observed

4. Pharmacopoeial Standards:

- DAC monograph requirements (where applicable)

Assay (99.0-101.0%):

- Based on analytical method variability (precision RSD ~1%)
- Consistent with ICH Q6A recommendations for APIs

Impurities (Any single $\leq 0.10\%$, Total $\leq 0.50\%$):

- Based on ICH Q3A(R2) qualification thresholds
- Manufacturing data: Maximum total impurities observed = 0.35%
- Provides adequate control margin

Water Content ($\leq 0.50\%$):

- Based on manufacturing data (typical: 0.25-0.35%)
- Non-hygroscopic compound (water uptake minimal)

Residual Solvent – Ethanol (≤ 5000 ppm):

- Per ICH Q3C Class 3 solvent (low toxicity)
- Typical levels: <500 ppm (well below limit)

3.2.S.5 REFERENCE STANDARDS OR MATERIALS

Primary Reference Standard:

- **Source:** European Pharmacopoeia (EDQM)
- **Name:** Fampridin for system suitability CRS (Chemical Reference Substance)
- **Batch Number:** [Batch No.]
- **Purity:** $\geq 99.5\%$ (as stated in Certificate)
- **Storage:** 2-8°C, protected from light
- **Expiry/Retest Date:** [Date]

Certificate of Analysis: On file

Use:

- Identity testing (IR, HPLC retention time)
- Assay determination (quantification)
- System suitability testing

Qualification:

- Primary reference standard is certified by pharmacopoeia (no additional qualification required)

Secondary/Working Standards:

- Not used (primary standard used directly for routine testing)

3.2.S.6 CONTAINER CLOSURE SYSTEM

Primary Packaging:

Container Type: Polyethylene-lined steel drum

Liner Material:

- Food-grade High-Density Polyethylene (HDPE)
- Thickness: 200 µm
- Meets EU Regulation 10/2011 (plastic materials in contact with food)

Drum:

- Material: Steel (exterior)
- Capacity: 25 kg API per drum
- Closure: Tamper-evident lid with locking ring

Labeling:

- Product name: Fampridin
- Batch number
- Manufacturing date
- Net weight
- Storage conditions: "Store at 15-25°C, protect from light and moisture"
- "For manufacturing use only" (not for direct patient use)

Protection Provided:

- Moisture barrier (HDPE liner + sealed drum)
- Light protection (opaque drum)
- Physical protection (steel drum)

Compatibility:

- No interaction observed between Fampridin and packaging materials
- Stability studies conducted in proposed commercial packaging (25 kg drums)

Specifications for Packaging Materials:

- HDPE liner: Meets food-contact regulations, Certificate of Compliance on file
- Steel drums: Clean, dry, free from rust or contamination

3.2.S.7 STABILITY

3.2.S.7.1 STABILITY SUMMARY AND CONCLUSIONS

Stability Studies Conducted:

1. Long-Term Stability (25°C / 60% RH):

- **Duration:** 24 months (ongoing to 36 months)
- **Storage Conditions:** 25°C ± 2°C / 60% RH ± 5% RH
- **Packaging:** Polyethylene-lined steel drums (25 kg)
- **Number of Batches:** 3 batches (pilot and commercial scale)

2. Accelerated Stability (40°C / 75% RH):

- **Duration:** 6 months
- **Storage Conditions:** 40°C ± 2°C / 75% RH ± 5% RH
- **Packaging:** Polyethylene-lined steel drums
- **Number of Batches:** 3 batches

3. Intermediate Stability (30°C / 65% RH):

- **Duration:** 12 months
- **Storage Conditions:** 30°C ± 2°C / 65% RH ± 5% RH
- **Packaging:** Polyethylene-lined steel drums
- **Number of Batches:** 3 batches

3.2.S.7.2 POST-APPROVAL STABILITY PROTOCOL AND STABILITY COMMITMENT

Ongoing Stability Program:

- At least **one batch per year** will be placed on long-term stability (25°C/60% RH)
- Stability data will be reviewed annually as part of Annual Product Quality Review (APQR)

Post-Approval Commitment:

- Continue long-term stability studies to 36 months (if not yet completed at time of approval)
- Report any significant changes or out-of-specification results to PMDA

3.2.S.7.3 STABILITY DATA

Summary of Stability Results:

LONG-TERM STABILITY (25°C / 60% RH) – Batch 2024-01-001

Time Point	Appearance	Assay (%)	Total Impurities (%)	Water (%)	Conclusion
Initial (T=0)	Conforms	99.8	0.32	0.25	Pass
3 months	Conforms	99.7	0.33	0.26	Pass
6 months	Conforms	99.7	0.34	0.28	Pass
9 months	Conforms	99.6	0.35	0.29	Pass
12 months	Conforms	99.6	0.36	0.30	Pass
18 months	Conforms	99.5	0.38	0.32	Pass
24 months	Conforms	99.5	0.40	0.33	Pass

Trend: Slight increase in impurities and water content over time, but all within specification limits.

ACCELERATED STABILITY (40°C / 75% RH) – Batch 2024-01-001

Time Point	Appearance	Assay (%)	Total Impurities (%)	Water (%)	Conclusion
Initial (T=0)	Conforms	99.8	0.32	0.25	Pass
1 month	Conforms	99.7	0.34	0.28	Pass
3 months	Conforms	99.6	0.38	0.32	Pass
6 months	Conforms	99.5	0.42	0.36	Pass

Conclusion: No significant change observed under accelerated conditions. All parameters within specification.

3.2.S.7.4 PROPOSED RETEST DATE / EXPIRY DATE

Based on stability data:

Proposed Retest Period: 36 months

Storage Conditions: Store at 15-25°C, protect from light and moisture

Justification:

- Long-term stability data (24 months) shows no significant degradation
- Accelerated stability data (6 months) shows product is stable
- Extrapolation to 36 months is justified per ICH Q1E guidelines
- Ongoing stability studies will continue to support retest period

Label Storage Statement: "Store at 15-25°C. Protect from light and moisture. Retest date: [36 months from manufacture date]"

APPENDICES

APPENDIX 1: ANALYTICAL PROCEDURES (ATTACHED SEPARATELY)

- HPLC Method for Assay and Impurities (HPLC-001)
- GC Method for Residual Solvents (GC-001)
- Particle Size Method (Laser Diffraction)

APPENDIX 2: VALIDATION REPORTS (ATTACHED SEPARATELY)

- AMV-HPLC-001: Analytical Method Validation Report – HPLC

- AMV-GC-001: Analytical Method Validation Report – GC

APPENDIX 3: BATCH MANUFACTURING RECORDS (SAMPLE) (ATTACHED SEPARATELY)

- Executed Batch Record for Batch 2025-01-001

APPENDIX 4: CERTIFICATES OF ANALYSIS (ATTACHED SEPARATELY)

- CoA for batches 2025-01-001 through 2025-05-001

APPENDIX 5: STABILITY DATA (DETAILED) (ATTACHED SEPARATELY)

- Full stability data tables and chromatograms

APPENDIX 6: SPECIFICATIONS FOR RAW MATERIALS (ATTACHED SEPARATELY)

- Starting Material A Specification
- Starting Material B Specification
- Ethanol Specification

APPENDIX 7: SUPPLIER QUALIFICATION REPORTS (ATTACHED SEPARATELY)

APPENDIX 8: EQUIPMENT QUALIFICATION REPORTS (SUMMARY) (ATTACHED SEPARATELY)

- IQ/OQ/PQ Summary for Reactors, Dryer, Mill, HPLC

APPENDIX 9: PROCESS VALIDATION REPORT (ATTACHED SEPARATELY)

- PVR-FAM-001: Process Validation Report

APPENDIX 10: IMPURITY QUALIFICATION (TOXICOLOGICAL ASSESSMENT) (ATTACHED SEPARATELY)

- Toxicological assessment for Impurity B

ABBREVIATIONS

- API: Active Pharmaceutical Ingredient
- CAPA: Corrective and Preventive Action
- CFU: Colony Forming Units
- CoA: Certificate of Analysis
- CPP: Critical Process Parameter
- CQA: Critical Quality Attribute
- CRS: Chemical Reference Substance
- CTD: Common Technical Document
- DAC: Deutscher Arzneimittel-Codex
- DMF: Drug Master File
- DoE: Design of Experiments
- EDQM: European Directorate for the Quality of Medicines
- GC: Gas Chromatography
- GMP: Good Manufacturing Practice
- HDPE: High-Density Polyethylene
- HPLC: High-Performance Liquid Chromatography
- ICH: International Council for Harmonisation
- INN: International Nonproprietary Name
- IPC: In-Process Control
- IQ/OQ/PQ: Installation/Operational/Performance Qualification
- IR: Infrared Spectroscopy
- JP: Japanese Pharmacopoeia
- LOD: Loss on Drying / Limit of Detection

- LOQ: Limit of Quantitation
- MS: Mass Spectrometry
- NMR: Nuclear Magnetic Resonance
- Ph. Eur.: European Pharmacopoeia
- PMDA: Pharmaceuticals and Medical Devices Agency (Japan)
- PXRD: Powder X-Ray Diffraction
- RH: Relative Humidity
- RSD: Relative Standard Deviation
- SEM: Scanning Electron Microscopy
- SOP: Standard Operating Procedure
- USP: United States Pharmacopeia

DOCUMENT APPROVAL

Prepared by:

[Name, Title: Regulatory Affairs Manager]

Signature: _____

Date: 04.12.2025

Reviewed by:

[Name, Title: QA Manager]

Signature: _____

Date: 04.12.2025

Reviewed by:

[Name, Title: Production Manager]

Signature: _____

Date: 04.12.2025

Approved by:

[Name, Title: Management Representative]

Signature: _____

Date: 04.12.2025

END OF DMF MODULE 3.2.S

Submission to PMDA: [Date To Be Determined]

Supporting Documents: See Appendices 1-10 (attached separately)